

Treadmill 7.0T

SPIRIT
powered by *dyaco*



The 7.0T treadmill is engineered to deliver superior rehabilitation and patient results. With a comprehensive suite of features, users can enhance their therapy sessions like never before.

This treadmill provides exceptional flexibility, allowing you to tailor every session to meet each client's unique needs:

- Full-length, adjustable height handrails for added stability
- -10% - 15% decline/incline range
- Low starting speed and step-up block for easy on/off access
- 3.0 HP AC motor with bidirectional walking direction capabilities
- Current Leakage Protection with EMC Filter

The 7.0T is designed to withstand intense therapeutic sessions, making it the ideal choice for any commercial therapy setting.



Intuitive, easy-to-read display



Fully adjustable parallel bars support a wide variety of users



Front and rear motors allow for true up/down walking

Specifications

Incline	-10% to 15%
Bidirectional Speed	Forward 0.1 - 10 mph and reverse 0.1 - 3 mph in 0.1 increments
Motor	3.0HP AC, IEC Class S6
Belt and Deck	22" x 60" self-lubricating deck
Power	120 volts, 60 Hz, dedicated 20-amp circuit (NEMA 5-20P standard power plug/receptacle type)
Max User Weight	440 lb / 200 kg
Dimensions	94.8" (83.7" w/o rear step) x 48.7" x 57.5"
Product Weight	496 lb
Current Leakage Protection	Yes, with EMC Filter
Compliance	CAN/CAS-C22.2 No. 60601-1:14, ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014, EN 60601-1-2:2015 IEC 60601-1:2015+A1:2012, EN 60601-1:2006+A1:2013+A12:2014 EN ISO 20957:2013, MDD 93/42/EEC Class IIa, NB: 0123
Warranty - Commercial (Dues paying) Rehab / Light Commercial	5 years parts, 1 years labor



NEMA 5 - 20P