

**MODEL A9069**

Shown with Optional 69 Gas Spring Backrest

Model A9069
6-Leg Dual Function Split Leg Table

- Extremely versatile and functional split leg design.
 - Each padded leg rest can be adjusted with knee fully extended and locked into desired angle up to 45° with spring-loaded lock-in lever.
 - Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.
- (6) legs and plywood top for added strength.
- Laminate shelf with (2) long stretchers.
- 30" x 78" upholstered top with 2" high-density foam.
- High-pressure laminate legs, apron and H-brace.
- New enhanced leg hardware system.
- Legs have protective black vinyl cuff caps.
- Natural Oak or Wild Cherry laminate.
- Choice of 9 Pro-Form™ vinyl upholstery colors.
- Full subtop for increased strength.
- Available with optional 69 Gas Spring Backrest. (400 lbs. weight capacity)



MODEL	LENGTH	WIDTH	HEIGHT
A9069	78"	30"	32"