



MODEL A9071

Shown with Optional 69 Gas Spring Backrest,
(3) A946 Lock/Bolts

**Model A9071
Dual Function Split Leg Lift Table**



- Extremely versatile and functional split leg design.
 - Each padded leg rest can be adjusted with knee fully extended and locked into desired angle up to 45° with spring-loaded lock-in lever.
 - Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.
- (6) Hinged doors with interior shelf offer concealed supply storage.
- Plywood top for added strength.
- 30" x 78" upholstered top with 2" high-density foam.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- Choice of 9 Pro-Form™ vinyl upholstery colors.
- Full subtop for increased strength.
- Available Options: 69 Gas Spring Backrest (400 lbs. weight capacity). 99L Left Hand Table Model, A946 Lock/Bolt for Double Doors.



MODEL	LENGTH	WIDTH	HEIGHT
A9071	78"	30"	31"