

# AFTR DC™ (Dorsal Closure)



## Ankle & Foot - Trauma Rehabilitation

1

Post-Surgical or  
Injury Edema Control

2

Compression  
Inside Walking Boot

3

Stability During  
Rehabilitation

The AFTR DC brace combines variable compression and stabilization to reduce pain, increase ROM, and improve recovery time.

### FEATURES

- *Variable Compression for Edema Control*
- *Dorsal Closure for Easy Post-Surgical Application*
- *Non-Elastic Stirrup Straps for Stability*
- *Silicone Gel Pads for Extra Edema Control Around the Malleoli*
- *Trimmable in Length for Perfect Fit*



<i>Product</i>	<i>Product #</i>
AFTR DC	5386*

**Note:** To indicate size, **replace \*** with number from price list.

## Edema Control

- **After Traumatic Injury or Surgery:** Dorsal closure facilitates donning after injury or surgery to immediately control edema and expedite return to normal range of motion.
- **Use with Walking Boot:** Allows foot to maintain necessary compression while comfortably fitting inside a walking boot.
- **Adjustable Compression:** Variable closure adapts to provide consistent compression during changes in foot volume.

## Ankle Stability

- **Stirrup Straps:** Non-elastic stirrup straps provide medial and lateral stability -- perfect step-down treatment after a walking boot.
- **Thin and Strong:** Durable material for daily use and thin enough to fit inside most shoes.

## Bio Skin Material

- **Comfortable:** Highly-compressive, tri-laminate material that is soft against the skin. Neoprene free and latex free.
- **Velcro Compatible:** The outside of the brace is Velcro-compatible to allow attachment of stirrup straps and WEIL accessories.
- **Breathable:** Heat and moisture vapor pass through the material.
- **Trimmable Fit:** Trim with scissors. Material will not fray.

## WEIL Digit Strapping Options

- **HCST™ (Hallux Control Strap):** For splinting and strengthening the hallux after bunionectomy.
- **WEIL Osteotomy Strap:** Prevents dorsal drift and floating toes.

## HCPCS

- L1902



Adjustable compression with easy dorsal closure



Stirrup straps offer ankle stability



Fits easily into most shoes



Weil Straps