WHIRLPOOL TABLES FEATURING STARBOARD ST TOPS



MODEL AT820 WHIRLPOOL TABLE FEATURING

MODEL AT825



STARBOARD® ST TOP

- Natural finish hardwood table with StarBoard® ST top and steps.
- StarBoard[®] ST is a high-density polyethylene sheet which is waterproof and scratch-resistant and features a textured surface.
- Standard color is Sanshade.
- Water-resistant urethane finish.
- Leg levelers for uneven floors.
- 47" wide x 32" deep available for 22", 24" and 26" whirlpools.
- Please specify manufacturer, model number and size and height of tank when ordering. Custom sizes available.

We add approximately $1 \frac{1}{4}$ " to the width of the tank as a clearance factor.

MODEL	LENGTH	WIDTH	HEIGHT
AT820 (not shown)	47″	32″	Variable

O/A footprint including steps is 56"L x 42"W

PROTEAM BY HAUSMANN • 130 Union St, Northvale NJ 07647 USA • 1-888-HAUSMANN • www.ProteamTables.com

MODEL AT825 300 WHIRLPOOL TABLE FEATURING STARBOARD® ST TOP WITH CABINET

- Same as AT820 but includes 2 shelf cabinet on backside of table.
- Shelf dimensions are 35"W x 12"H x 9"D.

MODEL	LENGTH	WIDTH	HEIGHT	
AT825	47″	32″	Variable	_
	1	// /0//\//		_

O/A footprint including steps is 56"L x 42"W

