

BARIATRIC ELECTRIC STAND-IN TABLE W/ PATIENT LIFT



Model 6177

- Motorized height adjustment of top from 42 ½" to 52 ½".
- Motorized patient lift can elevate a 500 lb. patient from sitting to standing position.
- Durable, white leather, high pressure laminate top 38" wide x 24.75" deep.
- Cutout size is 24" wide x 10" deep.
- Hand control operates electric table-top height and motorized patient lift.
- Includes lift support harness for standard and bariatric patients.
- Padded knee, chest and back supports.
- Heavy gauge steel frame in cream-powder coated finish.
- Frame base is 42" wide x 43" deep.
- Four (4) locking 3" diameter heavy-duty casters
- Wheelchair accessible.
- 110V, 60 Hz, 1.0 amp.



MODEL	WIDTH	DEPTH	HEIGHT
6177	43"	42"	42 ½" - 52 ½"