




BioSkin™ Extended Sizing Chart


Extended sizing is subject to a price increase and is nonrefundable. If interested in custom bracing, please call for more information.

Foot & Ankle

		By shoe size			
		TriLok™			
		Size	*PN	Men's	Women's
		S	1	5 – 6.5	6 – 7.5
		M	2	7 – 9.5	8 – 11
		L	3	10 – 12.5	11 – 13.5
		XL	4	13 – 16	14 – 16
		HCS - Hallux Control Strap			
		S	1	5 – 10	6 – 11
		L	3	10.5 – 16	11.5+ – 16

		Circumference from heel to top of the foot – back to heel			
		Size	*PN	Inches	Centimeters
		AFTR Products			
		XS - S	1	11 – 13	27.9 – 33
		M - L	3	12.5 – 14	31.8 – 35.6
		XL - 2XL	5	13.5 – 15	34.3 – 38.1
		Standard & Visco			
		XS	0	10 – 11	25 – 28
		S	1	11 – 12	28 – 31
		M	2	12 – 13	31 – 33
		L	3	13 – 14	33 – 36
		XL	4	14 – 15	36 – 38
		2XL	5	15 – 16	38 – 41
Extended Sizes		3XL	6	16 – 17	41 – 43

Knee

		Knee Products			
		Circumference with knee fully extended			
		Size	*PN	Inches	Centimeters
				mid – patella	7" above
				mid – patella	18cm above
		XS	0	12 – 13	16 – 18
		S	1	13 – 14	18 – 20
		M	2	14 – 16	20 – 22
		L	3	16 – 17	22 – 25
		XL	4	17 – 19	25 – 28
		2XL	5	19 – 21	28 – 32
		3XL	6	21 – 23	32 – 36
		4XL	7	23 – 25	36 – 40
		5XL	8	25 – 27	40 – 44
Extended Sizes				54 – 58	81 – 91
				58 – 64	91 – 102
				64 – 69	102 – 112

Indicate Sizing When Ordering

Product#	Size#	To indicate size, replace * with corresponding *PN (Part Number) from the sizing chart.
XXXXX	*	

Four Easy Ways to Order



ON-LINE

Pro.BioSkin.com

Easy, Fast, Secure, 24/7

- Place Orders
- Track Shipments
- View Account Balance
- Pay Bills



EMAIL

Orders@BioSkin.com

Email orders 24/7

Use Subject Line: "Bio Skin Order"



PHONE

1-800-541-2455

Our Customer Care team is ready to help.

Monday - Friday, 7am - 5pm Pacific time.



FAX

1-541-482-2119 (Fax orders 24/7)



Extended Sizing Chart

Wrist


Circumference of extended wrist				
Size	*PN	Inches	Centimeters	
XS - S	1	4 - 6	10 - 15	
M - L	3	6 - 8	15 - 20	
XL - 2XL	5	8 - 10	20 - 25	

Thumb Spica				
Length of thumb				
Size	*PN	Inches	Centimeters	
XS - M	2	3.25 - 4.25	8 - 11	
L - 2XL	5	4.25 - 5.25	11 - 14	

Thigh

Mid - thigh circumference				
Size	*PN	Inches	Centimeters	
XS	0	16 - 18	41 - 46	
S	1	18 - 20	46 - 51	
M	2	20 - 22	51 - 56	
L	3	22 - 25	56 - 64	
XL	4	25 - 28	64 - 71	
2XL	5	28 - 32	71 - 81	
Extended Sizes	3XL	6	32 - 36	81 - 91

Spine & Shorts

	Waist measurement			
	Size	*PN	Inches	Centimeters
	Baja™			
	S	1	26 – 30	66 – 76
	M	2	30 – 36	76 – 91
	L	3	36 – 44	91 – 112
	XL	4	44 – 54	112 – 137
	UNIV	4U	26 – 54	66 – 137
	Back Flex & Back Plus			
	1	1	28 – 34	71 – 86
	3	3	34 – 44	86 – 112
	5	5	44 – 56	112 – 142
	Back Skin			
	XS	0	24 – 26	61 – 66
	S	1	26 – 30	66 – 76
	M	2	30 – 34	76 – 86
	L	3	34 – 38	86 – 97
	XL	4	38 – 42	97 – 107
	2XL	5	42 – 46	107 – 117
	Extended Sizes	3XL	6	46 - 50
4XL		7	50 - 54	127 - 137
5XL		8	54 - 58	137 - 147
	Compression Shorts			
	S	1	28 – 32	71 – 81
	M	2	32 – 36	81 – 92
	L	3	36 – 40	92 – 102
	XL	4	40 – 44	102 – 112
	2XL	5	44 – 48	112 – 122
	Extended Sizes	3XL	6	48 - 52
4XI		7	52 - 56	132 - 142

Elbow

Circumference of largest part of forearm				
Size	*PN	Inches	Centimeters	
XS	0	8 - 9	20 - 23	
S	1	9 - 10	23 - 26	
M	2	10 - 11	26 - 28	
L	3	11 - 12	28 - 31	
XL	4	12 - 13	31 - 33	
2XL	5	13 - 14	33 - 36	
Extended Sizes	3XL	6	14 - 15	36 - 38
	4XL	7	15 - 16	38 - 41
	5XL	8	16 - 17	41 - 43

Calf

Largest part of calf circumference				
Size	*PN	Inches	Centimeters	
S	1	11 - 12	28 - 30	
M	2	12 - 14	30 - 35	
L	3	14 - 16	35 - 41	
XL	4	16 - 18	41 - 46	
2XL	5	18 - 20	46 - 51	
Extended Sizes	3XL	6	20 - 22	51 - 56