



ACTIVE INNOVATIVE  
 CONFIDENT BE BETTER  
 PERFORMANCE DYNAMIC COMFORT



**NEW** SlimLok™ Adjustable Hinge


Featured in the Q Lok ROM and CrossFire ROM

- **Thinner**  
One of the thinnest adjustable hinges on the market. Only 1/2 inch thick with multiple flexion and extension stops.
- **Easy Flexion/Extension Adjustment**  
Adjusting the hinge angle is easy with the extension slider and flexion dials. No loose wedges or specialty tools needed. Easily adjustable with a penny.
- **Flared Hinge Uprights**  
Metal hinges are flared to conform to a leg's natural contour for improved comfort and stability.

**NEW** ActiveSeam™ The Stitch Matters®

Featured in the Q Lok ROM and CrossFire ROM

- With no overlapping material, the new ActiveSeam is almost undetectable against the skin.
- **30% Stronger than a Flatlock stitch**
- **Flatter, smoother seam.**
- **Twice as much stretch as a Flatlock stitch**



**NEW** Cirrus™ Material

Featured in the Q Lok ROM and CrossFire ROM

- **Ultra-Thin**  
Breathable, single-layer material with high-performance stretch.
- **Enhanced Comfort**  
Perfectly conforms to the knee providing maximum comfort and freedom of movement.
- **No Bunching**  
Eliminates popliteal bunching and irritation.

Products

	Part Number
Q Lok ROM	6166*
Crossfire ROM	6180*
Visco Knee Skin w/ Conforma Hinge - Stratus Material	6137*
Patella Stabilizer w/ Conforma Hinge - Stratus Material	6132*
Patella Stabilizer w/ Conforma Hinge - Ultima Material	4132*
Q Brace w/ Conforma Hinge - Ultima Material	4136*
Q Brace w/ Conforma Hinge - Front Closure	5138
Premium J w/ Conforma Hinge - Stratus Material	6138*
Gladiator 15° ROM Hinge - Ultima Material - Pull-on	4134*
Gladiator 15° ROM Hinge - Stratus Material - Pull-on	6136*
Gladiator 15° ROM Hinge - Front Closure	5136*
Hinged Knee Skin - Ultima Material	4130*
Hinged Knee Skin - Stratus Material	6130*
Hinged Knee Skin - Front Closure - 2SL Material	5133*
Q Lok Patella Traction Brace	6164*
Q Lok Patella Traction Brace - Front Closure	5164*

Note: To indicate size, replace \* with corresponding number from the sizing chart.

Sizing

		Circumference with knee fully extended			
Size	*PN	Inches		Centimeters	
		mid - patella	7" above	mid - patella	18cm above
XS	0	12 - 13	16 - 18	31 - 33	41 - 46
S	1	13 - 14	18 - 20	33 - 36	46 - 51
M	2	14 - 16	20 - 22	36 - 41	51 - 56
L	3	16 - 17	22 - 25	41 - 43	56 - 64
XL	4	17 - 19	25 - 28	43 - 46	64 - 71
XXL	5	19 - 21	28 - 32	46 - 54	71 - 81



240 East Hersey St., Suite 2 Ashland, OR 97520-5202  
 t: 800-541-2455 f: 541-482-2119

CropperMed.com



Patent #: 5,735,807  
 ©2015 Cropper Medical, Inc.  
 03/2015 PN: 12475A



**BioSkin™**

Be Better™

## Our Most Comfortable Braces Ever!

The CrossFire ROM™ and the Q Lok ROM™ feature Bio Skin's new, ultra-thin Cirrus™ material and our new ActiveSeam™ stitch. Both products are constructed with a new cut that prevents migration and eliminates bunching behind the knee, making these our most comfortable braces ever.

NEW



### CrossFire ROM™ Knee Ligament Protection Featuring the New SlimLok™ Adjustable Hinge

The CrossFire ROM is designed to protect knees with ligament tears or sprains. The rigid uprights are slightly flared to contour to the leg's natural shape and provide stability so patients can move and be active with confidence. The thin hinge is easy to adjust for flexion and extension control during patient rehabilitation.

#### FEATURES

- Length: 13"
- Fits Left or Right
- Easily Adjustable Flexion / Extension Settings
- Available in Pull-On or Front-Closure Design
- ActiveSeam Stitch

#### INDICATIONS

- Ligament Sprains / Tears
- Meniscus Tears
- Knee Instability
- HCPCS: L1832 / L1833

NEW



### Q Lok ROM™ Patella Traction System Featuring the New SlimLok™ Adjustable Hinge

The Q Lok ROM is designed to diminish the effects and aid in the rehabilitation of patellofemoral pain. The Dynamic Patella Traction System exerts a medial pull on the patella to improve tracking in the trochlear notch and increase surface contact area in the patellofemoral joint.

#### FEATURES

- Length: 13"
- Addresses Glide and Tilt
- Fits Left or Right
- Easily Adjustable Flexion / Extension Settings
- Available in Pull-On or Front-Closure Design
- ActiveSeam Stitch

#### INDICATIONS

- Anterior Knee Pain
- Maltracking Patella
- Patellofemoral Osteoarthritis
- Chondromalacia
- HCPCS: L1832 / L1833



### Visco Knee Skin w/ Conforma Hinge

Functional compression for swelling and pain.

#### FEATURES

- Silicone Gel Ring
- Finger Loops for Easy Application
- Closed Patella

#### INDICATIONS

- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis
- HCPCS: L1820



### Q Brace™ w/ Conforma Hinge

Adjustable support for patellofemoral pain.

#### FEATURES

- Fits Left or Right
- Available in Pull-On or Front-Closure Design

#### INDICATIONS

- Anterior Knee Pain
- Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain
- HCPCS: L1820



### Patella Stabilizer w/ Conforma Hinge

Compression and support for mild to moderate knee pain.

#### FEATURES

- Foam Ring Around Patella for Stability

#### INDICATIONS

- Patella Stabilization
- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis
- HCPCS: L1820



### Premium J™ w/ Conforma Hinge

Simple support for patellofemoral pain.

#### FEATURES

- Adjustable, Lateral Patellar Tension
- Closed Patella
- Specify Left or Right

#### INDICATIONS

- Anterior Knee Pain
- Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain
- HCPCS: L1820



### Gladiator Knee Ligament Protection

Maximum stability for damaged ligaments.

#### FEATURES

- Length: 15"
- Silicone Gel Patella Ring
- Easily Adjustable ROM Settings
- Open Patella
- Available in Pull-On or Front-Closure Design

#### INDICATIONS

- Ligament Tears or Sprains
- Knee Instability
- HCPCS: L1832 / L1833

NEW HINGE LIGHTER & THINNER



### Hinged Knee Skin Lightweight & Low-Profile

Moderate support for damaged ligaments.

#### FEATURES

- New Hinge is 33% Lighter & 40% Thinner
- Finger Loops for Easy Application
- Strapping System Minimizes Migration
- Available in Pull-On or Front-Closure Design

#### INDICATIONS

- Ligament Tears or Sprains
- Meniscus Tear
- HCPCS: L1820



### Q Lok™ Dynamic Patella Traction

A system that decreases pain by increasing surface contact area in the joint.

#### FEATURES

- Addresses Glide and Tilt
- Fits Left or Right
- Lightweight: 7 Ounces
- Available in Pull-On or Front-Closure Design

#### INDICATIONS

- Patellar Tilt and Glide
- Patellar Subluxation
- Post Surgical Lateral Release
- Patellofemoral Osteoarthritis
- HCPCS: L1820