

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

ILLUVIA LASER

STANDARD FEATURES:

- Unit Dimensions: 12.5" W x 10.75" L x 4.5" H
- Interface control panel is a color LCD touch screen
- Folding touch screen for ease of use
- Operation mode: CW-single-repeat
- Quick start or pre-set protocols
- Accessories: 1 Hard plastic rolling carry case, DC 15V adapter, 3 laser safety goggles, 1 set of PM5 therapy handpiece with 5 treating probes (30mm big massage ball, 30mm non contact, 30mm contact, 5mm accupuncture, 10mm contact)



GBOX-5W

Class 4 Laser according to EN/IEC 60825-1:2014

LASER CONFIGURATIONS:

GBOX-5W - Illuvia 5W (810 nm)

QTY: _____

GBOX-30W - Illuvia 30W (810 nm + 980 nm)

QTY: _____



*Illuvia Carry
Case*



*Illuvia Carry
Case (open)*