




KoreSway Kompact

Model Number: KoreSway Kompact-19



We're driven because you're driven. The KoreSway™ Kompact is designed for use by medical centers and clinics, medical and PT offices, assisted living facilities, recreation centers, health and fitness clubs, universities and professional sport teams.

A Study in Perfection. The Korebalance™ System uses the latest in virtual and interactive technology, offering innovative high-tech balance assessment and training. Backed by over two decades of medical experience, the Korebalance™ systems will provide you and your patients with the best possible documented neurosensory assessment and training.

SPECIAL FEATURES:

- 19" or 22" Touchscreen Computer
- Fixed Monitor Mount
- Fixed Safety Rails
- Force Plate Platform with 23" (58.4 cm) Diameter
- Includes Korebalance® testing/training software suite and games
- Optional Cognitive, Parkinson's Software Suite
- Patient Weight Capacity: 330 lbs (150 kg)
- Machine Weight Capacity: 112 lbs (50 kg)
- Dimensions: 27"x37"x67"
- Max Height: 73" (15.42 cm)
- Software Available in English, Spanish, French, German and Mandarin
- Frame Color: White
- Shipping Weight: 75 lbs (34 kg)
- Shipping Dimensions: 47" x 27" x 37"

Warranty: Frame: 5 year; Computer: 1 years; Parts: 1 year

In the interest of product development and improvement, specifications subject to change without notice.