

Midline

Positioning Wheelchair

POSITIONING
BRODA[®]




Full tilt and recline functionality provides superior positioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort





Comfort Tension Seating[®] increases support as well as sitting tolerance

Patient, Facility, & Caregiver Benefits


- Increases quality of life and engagement with surroundings
- Allows for easy maneuverability and safe patient adjustments
- Reduces risk of injury and falls with supportive stable positioning

Enhanced positioning with refined features

UPHOLSTERY OPTIONS



FRAME OPTIONS



OPTIONAL ACCESSORIES



Foot Box



O₂ Holder



Lateral Support



Half Tray



Tray



Padded Thigh Strap

SPECIFICATIONS

Degrees of Tilt	Up to 38°
Recline	Up to 90° (Trendelenburg capable)
Weight Capacity	up to 350 lbs. (114 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm)
Seat Depths	15.5", 17", 18.5", 20" (39cm, 43cm, 45cm, 58cm)
Floor-to-Seat Heights	Adjustable 17", 18", 19", 20", 21", 22" (43cm - 55cm)
Wheels	5" Front Casters: Total Lock (12.5cm) 5" Rear Casters: Directional or Total Lock (20cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating [®]
Standard Padding	Neckrest, shoulder bolsters, upper and lower side pads, seat and back pad, calf and sole pad