



HCI FITNESS

PhysioMill
Rehabilitation Treadmill

MSRP: \$5,499
Call Now!
800-398-9121
360-321-7090

CONTACT US



Ask About Our PhysioGait Packages!



The **PhysioMill** rehabilitation treadmill has a heavy duty 4 HP **Hi-Torque** DC Servo motor specifically designed to offer a true zero start with .10 mph increments for patients up to 500 lbs., with speeds up to 15 mph.

The extremely versatile **PhysioMill** rehabilitation treadmill has a **reverse belt** mode with 15 percent elevation allowing for downhill negative walking.

The 500 lb. weight capacity and the 20" x 63" walking belt allows the **PhysioMill** to be used in a variety of clinical settings such as Neurological, Orthopedic, Sports Medicine, Cardiac Rehab, Bariatrics and Pediatrics.

The **PhysioMill** rehabilitation treadmill comes standard with; parallel medical side rails, dual emergency stop buttons, hand

Click to Print PhysioMill Specifications

PhysioMill	
Technical Details	
Drive System	4 HP Continuous-Duty DC Servo Motor
Belt Direction	Foward and Reverse Direction
Running Surface	22" W x 63" L (56cm x 160 cm)
Belt	1" (2.5 cm) Thick Reversible Belt
Treadmill	Four-Ply ErgoDeck Shok Absorption System

Deck	
Treadmill Rollers	3.5" (8.9 cm), 22 lbs (10 kg) Steel Rollers
Speed Forward	0 - 12 mph (0 - 19.2 kph)*True Zero Start
Speed Backward	0 - 3 mph (0 - 4.8 kph) *True Zero Start
Speed Increments	0.1 Speed Adjustment Increments
Incline Range	0 - 15%
Max User Weight	500 lbs.
Deck Height	Only 6.25" (16.5 cm) Step-Up
Emergency Stop	Dual Emergency Stop Buttons & Safety Lanyard
Hand-Held Remote	Therapist Start-Stop
Domestic Power	120 VAC, 60Hz, 15 A
International Power	220 VAC, 50/60Hz, 15 A
Product Type	Commercial Grade
Treadmill Model #	TRM
HandRail Kit Model #	TRM-HR
Computer Details	
Display	High Definition 9" LCD Display

held remote start/stop/pause, orthopedic walking belt and contact heart rate.

The **PhysioMill** may be the last rehabilitation treadmill you will ever need to purchase. **“Over built and built to last!”**

Product Videos

[PhysioStep LXT - Linear](#)

MORE VIDEOS

Related Products



CMT 20-58

ALL PRODUCTS

Testimonials



Customer Reviews

Read product reviews and customer testimonials for HCI Fitness products!

Backlight	Easy to Read High Contrast Blue
USB Summary	Export Workout Data Including Peaks and Averages
Feedback	Speed, Time, Distance, Grade, Pulse, Calories, METS & Reverse
Heart Rate	Contact Heart Rate on Seat Handles & Built in Receiver Compatible with Polar and other Manufactures *Chest Belt Not Included
*Light Instiutional Warranty	
Frame	Lifetime
Mechanical Parts	5 Years
Electronic Parts	5 Years
Labor	1 Year
Shipping & Assembly	
Shipping	Freight
Assembly Required	Support Arms, Display, Side Rails
Unit Dimensions	
Compact Dimensions	60" L X 30" W X 40" H
Unit Weight	210 lbs.

*Excludes heavy use pay-for-membership facilities including health clubs, gyms, and YMCAs

PRODUCT BROCHURE

USER MANUAL

CONTACTSOURCES