PneuWeight

UNWEIGHTING SYSTEM

The PneuWeight Unweighting System allows patients to perform therapeutic activities in an upright, functional position without altering gait patterns in a safe and secure environment for both client and practitioner.

A TRULY FLEXIBLE SYSTEM

- Unique cantilevered design for unobstructed access.
- The largest range of unweighting in the industry: o to 300 lbs.
- 30" vertical oscillation offering Balance and Exercise modes, as well as a drop-stop.
- 360 degree rotation allows multi-planar exercise.
- Patented vest with adjustability for posture and gait.
- Adjustable height range from 7'11" to 8'8".
- Sturdy and simple to use control panel.
- Customized pricing available.
- Available in single and double stations.
- Available with Analog control panel.

Key Benefits

- Abilities patients
- **Bariatric** patients
- Orthopedic patients
- Strokes
- Back problems
- Gait training
- Cardiovascular
- Respiratory ailments
- Neurological patients
- Safety



Analog Control Panel Easy to read, easy to use analog control allows precise unweighting from o to 300 lbs. in 1 lb. increments. Features balance and exercise modes in familiar dial format.



All PneuWeight systems require a compressor to operate.



Patent 6672998

The system can be used over most exercise equipment, including but not limited to: treadmills, vibration platforms and balance equipment.

"I can't open my new clinic without the PneuWeight Unweighting System. It is one of the few capital purchases that I think is essential to my practice." Dr. Bill O'Grady, DPT, MA, OCS,

FAAOMPT, FAAPMa

"The PneuWeight system is the most effective program for dealing with treatment of back pain and prevention of injury. The results are consistent and reproducible between clinics." Steve Sodorff, PT

Unweighting Studies on the Web www.pneuthera.com/studies

Specifications



ACCESSORIES

One-PneuVest *

Your choice of style and size

Light weight, easy to use. Securely holds and lifts patients. Three styles:

Sizes: Vectors XSmall to Sport XLarge Abilities 3-in-1









One-PneuGait Assortment

One PneuGait Foot strap One PneuGait Drop Foot Strap One PneuGait Assist Foot Strap and Band







*Two-PneuVest for a Double Station

One-Cable Assembly

Extra safety with design that allows 360° movement.



One-Chain Extender Provides extra length on cable support bar while unweighting

One-Safety Drop

Limits the patient drop to 2" if they stumble



CONFIGURATIONS AVAILABLE:



Single Station

Includes One PneuVest



Double Station

Includes Two PneuVest

Specifications

	Single	Double
Height	7'11"- 8' 8"	7' 11" - 8' 8"
	242 - 265 cm	242 - 265 cm
Width	4' 2"	6' 4"
	127 cm	193 cm
Max	300 lb	300 lb*
Lifting	136 k	136 k
* Each station		

All PneuWeight systems require a compressor to operate.