

**NEW UNIVERSAL DESIGN:**  
*Fits either left or right leg*

# PREMIUM J™

## WITH CONFORMA HINGE

L1820

**A simple and comfortable solution for patellofemoral pain and instability**

- J-shaped silicon buttress forms a lateral wall to prevent subluxation and encourage proper alignment. Packaged with bilateral buttresses; simply remove one side to fit the left or right leg.
- Elastic, anterior straps control the patella to correct glide and tilt. May be positioned on the right or left.
- Also available without a hinge (AA4467).
- Available in sizes XS-5XL.



lightweight



dynamic stretch



breathable



hypoallergenic

PRODUCT

PART #

Premium J with Conforma Hinge

6138\*U

Premium J without hinge

6163\*U