



SMART



A Kismet Brand

stravasolutions.com

STRAVA SMART SPLINT BATH

A CLEAN, SMART, INNOVATIVE DESIGN

A state-of-the-art, first-of-its-kind splint bath precisely engineered to remove the guesswork and help prevent cross contamination

Automatic Dual Disinfection

The only splint bath designed with two built-in methods of disinfection to keep patients safe and meet regulatory compliance

Fast & Uniform Heating

Only 25 minutes to reach set water temperature with uniform and accurate heat distribution

Smart Technology

Preset and custom thermoplastic settings, auto on/off, low water sensor, alerts, water evacuation, and much more

kills
99.9%
of introduced
microbes

Meet Regulatory Compliance





- ✓ UV-C Germicidal Kill
- ✓ High-Heat Sanitation Cycle
- ✓ Thermoplastic Presets
- ✓ Two-Stage Filtration
- ✓ Auto On/Off
- ✓ Fast & Uniform Heating
- ✓ Audible/Visual Alerts
- ✓ Bi-Fold Stainless Steel Lid
- ✓ Low Water Sensor
- ✓ Rapid Water Evacuation

THERMOPLASTIC PRESETS

Recommended time and temperature presets




Choose pre loaded thermoplastic settings from top manufacturers or create your own from the LCD touch screen – only available on the Strava Splint Bath.



DIGITAL TOUCH INTERFACE

Traditional splint bath controls are antiquated and imprecise

A full menu of features allow clinicians control and customization from an easy to use touch display



Fully retractable bi-fold lid fits under most cabinets for efficient use of counter top space



STRAVA

info@stravasolutions.com
www.stravasolutions.com

**Ask Your Preferred Healthcare
Supplier For A Strava Splint Bath**



©2023 Strava Solutions



BASIC



A Kismet Brand

stravasolutions.com

INTRODUCING A NOT-SO-BASIC SPLINT BATH

NEW

The Strava Basic Splint Bath offers innovative features and modern looks with an economy price

Modern Looks


Made of a thermoset material, known for its insulation, flammability, resistance and durability. Features a stainless steel hinged bi-fold lid for easy access from three sides.

Auto On/Off

Automatically preheats the water to reach a set temperature at a set day and time. Ready when your clinic opens and energy efficient at night when not in use.

Heats in 25 minutes

Internal heat exchanger is immersed directly in the water basin, ensuring fast heating and uniform heat distribution.







- ✓ Auto On/Off
- ✓ Fast & Uniform Heating
- ✓ User Sets Time & Temperature
- ✓ Audible/Visual Alerts
- ✓ Low Water Sensor
- ✓ Rapid Water Evacuation
- ✓ Bi-Fold Stainless Steel Lid

DIGITAL TOUCH INTERFACE

Traditional splint bath controls are antiquated and imprecise

A full menu of features allow clinicians control Strava's features.



Fully retractable bi-fold lid fits under most cabinets for efficient use of counter top space





Economy Model


UVC Effective
Sterilization Model

STRAVA **BASIC** -OR- STRAVA **SMART**

Automatic On/Off
Heats in 25 Minutes
User Sets Time & Temp
Visual Timer
Temperature Display
Audible/Visual Alerts
Digital Touch Interface
Low Water Sensor
Bi-Fold Stainless Steel Lid
Rapid Water Evacuation

All the features of the Basic plus:

Water Circulation System
Dual Disinfection
UVC Germicidal Kill
High Heat Sanitation Cycle
Thermoplastic Presets
Custom Presets
Dual Filtration



kills
99.9%
of introduced
microbes



info@stravasolutions.com
www.stravasolutions.com

©2023 Strava Solutions