

Full length medical handlebars
(adult and pediatric heights)

5 HP permanent
magnet AC motor

Reverse speed
up to 3 MPH

3% decline to
12% incline



TR8000i - COMMERCIAL TREADMILL (MEDICAL)

CONSOLE FEATURES

| | |
|------------------|---|
| Display | Large 6.5" backlit LCD screen with scrolling alphanumeric messaging, program profile showing workout progress, and simultaneous display of all readouts Three additional LED displays with readouts for time, speed, and incline |
| Readouts | Program name, date and time, speed, incline, distance, calories, steps, heart rate and p |
| Programs | Two heart rate control programs, and 17 programs grouped by exercise motivation for weight loss, healthy living, and sports training, and Rockport Fitness Test |
| Quickstart | Yes |
| USB Ports | 2 USB ports to automatically save your exercise data and charge your wireless mobile device simultaneously |
| Bluetooth | Wirelessly connect to your supported iOS tablet for additional console functionality and control including graphical heart rate tracking, elevation change, and average heart r |
| Heart Rate | Grip pulse and heart rate telemetry |
| Media Holder | Book, magazine, and tablet computer holder |
| Intelli-Key™ | Button lighting system prompts users through program set-up and operations |
| Emergency Stop | Magnetic safety clip, when disengaged, immediately turns off power, bringing the treadmill to a gentle stop |
| Units of Measure | Standard and metric settings |
| Language Support | English |

MECHANICAL FEATURES

| | |
|-------------------|--|
| Speed | Forward: 0.1 to 12.0 MPH / 0.1 to 20.0 KPH Reverse: -0.1 to -3.0 MPH / -0.1 to -5.0 KPH |
| Incline | 3% decline to 12% incline |
| Running Surface | 22" W x 62" L / 56cm W x 157.5cm L |
| Belt | Low maintenance, multi-ply polyester belt with a PVC non-slip surface |
| Drive Motor | 5 HP permanent magnet AC |
| Incline Motor | 1,000 lbs. maximum lift |
| Rollers | 3" crown rollers keep belt centered |
| Deck Suspension | Independent shock absorbers |
| Deck | 1" melamine laminated reversible deck medium-density fiber board |
| Power Requirement | 110 - 120 VAC / 220 - 240 VAC (depending on country standard) |
| Dimensions | 80.5" L x 39.5" W x 61" H / 204.5cm L x 100.3cm W x 154.9cm H |
| Net Weight | 378 lbs. / 171.5 kg |
| Gross Weight | Box 1: 400 lbs. / 181.4 kg Box 2: 52 lbs / 23.6 kg |
| Max. User Weight | 500 lbs. / 226.8 kg |

| Warranty | Residential | Commercial <i>Up to 6 hours per day (TR7000i)</i> <i>Up to 10 hours per day (TR8000i)</i> |
|-----------------|--------------------|--|
| Frame: | Lifetime | 10 Years |
| Motor: | Lifetime | 5 Years |
| Parts: | 5 Years | 5 Years |
| Labor: | 2 Years | 1 Year |

LifeSpan

www.LifeSpanFitness.com

US Sales

Email: sales@LifeSpanFitness.com

Phone: 877.654.3837