


AutoPrism

New Light Therapy from Rich-Mar



- **More Power** - Over 600 mW of power at multiple wavelengths making it the most powerful, versatile, and economical light therapy device available.
- **Research-Based Effectiveness** - Proven red and infrared wavelengths to treat a wide variety of soft-tissue conditions.
- **Quick and Easy to Use** - Can be used as a handheld for quick 1J/cm² applications in 30 seconds or secure it in place for Hands-Free Light Therapy for longer protocols delivering 10 J/cm².
- **Affordable** – Separate, economical AutoPrism product as well as an AutoPrism plug-in accessory for AutoSound devices.

Rich-Mar Makes Waves in Therapy...

The new, innovative **AutoPrism Light Therapy** device provides both red and infrared wavelengths to treat a wide variety of soft-tissue conditions quickly and easily.

More Power – Treat larger areas in less time.

Rich-Mar announces another new innovation in therapy - Automatic Hands-Free Light Therapy, the AutoPrism. This unique, Automatic Light Therapy device provides over 600 mW of power at therapeutic wavelengths that will allow you to treat a wide variety of soft-tissue conditions.



Easy to Use - Handheld or Strap-on Application

It can be used as a handheld piece to deliver quick treatments of approximately 1 J/cm² in 30 seconds. It can also be used as an automatic, Hands-Free, strap-on light therapy applicator to deliver longer, research-based 4 J/cm². You can even use it for longer time protocols to deliver up to 10 J/cm². The unique Hands-Free use of the AutoPrism applicator allows you to use ALL of the various protocols to apply more energy to the patient without taking up more of your time.



Economical Separate Light Therapy Modality

Another great feature of the AutoPrism is that it can be purchased as a stand-alone modality. Therefore you do not have to purchase a new stimulator or combination device in order to add a Light Therapy device to your practice. Of course you will also be able to get it as an attachment for any of the AutoSound Hands-Free Ultrasounds and/or Combo devices.

For more information visit us on the web at www.richmarweb.com or call us at 800.762.4665.

Authorized representative:

