

Comfort

Rehab Wheelchair

REHAB
BRODA[®]



Shown with optional head support, Matrix Elite back support, Matrix Vi seat cushion, and swing-away lower leg supports

Patient, Facility, & Caregiver Benefits

- Pressure redistribution that aids in maintaining skin integrity
- Durable frame to withstand heavy workload healthcare environments
- Adaptable for wide range of patients and care needs



Front pivot seat tilt maintains user comfort with smooth, effortless tilt



Posterior tilt optimizes pressure redistribution



Solid seat pan offers support and is compatible with a variety of seat cushion options



Range of simple adjustments accommodates a variety of needs

Matrix[®] seating options available to ship preconfigured on-chair.

matrix[®] SEATING SERIES

Elan Headrest



Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.

Elite TR Back



Shell - Contoured aluminum with cutouts for lifting and breathability. Available in standard (3") OR deep (6") contours.

Foam - HR (high resilient) foam enclosed in cover.

Libra Cushion






Foam - HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.

Simplified Durability for Rehabilitation Care

SPECIFICATIONS

<i>Degrees of Tilt</i>	Up to 45°
<i>Weight Capacity</i>	Up to 250 lbs. (114 kg)
<i>Angle Adjustable Back Canes</i>	5°, 10°, 15°, 20°
<i>Seat Widths</i>	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
<i>Seat Depths</i>	Adjustable: 16", 17", 18", 19", 20" (40cm, 43cm, 45cm, 48cm, 50cm)
<i>Seat to Surface Heights</i>	Adjustable: 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19" (35cm, 38cm, 39cm, 42cm, 43cm, 44cm, 45cm, 48cm)
<i>Wheels</i>	5", 6" or 8" Front Caster Wheels with fork (13cm, 15cm, 20cm) 20", 22" or 24" Rear Mag Wheels (50cm, 55cm, 60cm)
<i>Lifetime Warranty Frame</i>	10 years on 16 Gauge tubular steel frame
<i>Government Coding (US & CA)</i>	HCPCS: E1161 ADP (ON): WMA5B-025

FRAME OPTIONS



Additional accessories available: lower leg supports, trays, and positioning belts